Modelo 146 M, biológico | Model 146 M, biological

Marca|Trademark: Zuzi

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The Zuzi microscope model 146M is industrial type microscope with episcopic and diascopic illumination that uses an infinity corrected optics assuring an excellent quality of the images obtained and presenting a wide range of applications in metallurgy, mineralogy and materials science for studying minerals, PCBs, films, liquid crystals, micro-electronics, etc. A very robust instrument with a T design that guarantees a good stability and the absence of vibrations that could disturb observation; without losing ergonomics, with knobs at a lower position and a Siedentopf-type head to provide the user a comfortable and tension-free position. The microscope uses infinity corrected plan achromatic objectives with larger working distances than usual objectives thus being specially appropriated when observing thicker samples.

Observations can be performed under reflected or transmitted light with a complete light intensity control by the user thanks to the illumination Köhler system:

■ Reflected light system: the halogen lamp is located at a compartment and presents centering knobs. The light intensity can be adjusted by the field and aperture diaphragms and a rotary filter-holder that includes green, yellow and blue filters.

To perform observations under polarized light the microscope is equipped with a polarizer placed at the filter holder and an analyzer rotary in 360° that can removed from the light path.

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■ Transmitted light system: for studying transparent or translucent materials the base of the microscope incorporates an halogen lamp. With height-adjustable Abbe condenser (N.A. 1.25), field and aperture diaphragms, and mobile filter holder that includes green, blue and yellow filters. The trinocular tube is compatible with Zuzi eyepiece video cameras allowing capturing and recording the images obtained, thus increasing the equipment applications and turning it into a complete and versatile microscope.

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